

Fractional CO2 laser



Advantages

RF driving tube, no need to change. It emits very fine and well arranged laser which can penetrate to corium layer directly.

Best treatment result which can replace surgery totally. No wound, quick time convalescence. Wide range of application, effective treatment.

8.4 inch real LCD color touch screen, most humanized software control
7-jointed articular light conduct, accurate treatment. No pain, no cut, no side-effects! No rebound! Apparent long-last result!

Specifications	Fractional CO2 laser	
LCD screen	8.4 inch LCD real color touch screen	
Light conducting system	7 jointed arms to conduct light	
Wavelength	10.6µm	
Pulsed Radiofrequency	0.530W	
Power	30W	
MODE	Pulse	Scan
Laser energy	1-30W	5-45mj
Step value	1W	5mj
Scan pattern size	0.1x0.1mm—20x20mm	
Scan Mode	free scan; Sequence scan; maximum distance scan	
Dot quantity	400 dots at most	
Pulse duration	0.1-10ms	
Spot distance	0.1-2mm	
Pulse energy	1mj to 100mj is optional for each dot.	
Cooling system	Water and wind cooling system	
Laser apparatus	Sealed off laser device stimulated by direct current	
Condenser focus	f=50mm	
Spot size	≥0.2mm at the focus	
Condenser focus	f=50mm	
Radiation time/ Interval time	0.01 1 s	
Aiming beam	635nm infrared ray	
Power supply	110/220V + 10%50Hz/60HZ 6A	
Continuous working period	24 hours	
Language	English, Spanish, French, Portuguese, German, Russian	
Price	Click here!	
Please ask!		